

Small and practical **C-Arm**



CONTACT US

INTERMEDICAL S.r.l.
IMD Group

24050 Grassobbio (Bergamo) Italy
Via E. Fermi, 26
ph. 0039 035 6594811
fax 0039 035 6594899
e-mail info@inter-med.it
www.inter-med.it

Data may be altered with no prior warning.

doc.12983 - 04.11.2019



New Radius
Single




IMPROVING PERFORMANCES

Membrane keyboard with alphanumeric touch screen LCD display for all the parameters and error messages.
Keyboard can be rotated of $\pm 60^\circ$.

*Tastiera a membrana con schermo LCD touch-screen per la gestione di tutti i parametri operativi e i messaggi di errore.
La tastiera può essere ruotata di $\pm 60^\circ$.*

!!! WARNING !!!
This X-ray unit may be dangerous to patient and operator unless safe exposure factors, operating instructions and maintenance schedules are observed



MONITOR SUPPORT MOVEMENTS.

Movimentazioni braccio monitor.



A small and practical C-Arm with one orientable 24" LCD touchscreen monitor, assembled directly on the mobile equipment (no trolley).

Stativo mobile ad arco a 'C', piccolo e pratico, con un monitor LCD da 24", touchscreen, orientabile, posizionato direttamente a bordo macchina (senza carrello).

MODELS AVAILABLE:

with 9" Image Intensifier and fixed or rotating anode monoblock.

MODELLI DISPONIBILI:

Con amplificatore di brillantezza da 9" e monoblocco con anodo fisso o rotante.

NEW DIM HD SOFTWARE (EYES):

Hard Disk SSD 128 GB (55.000 images) Measure Software included DICOM available as options

NUOVO SOFTWARE DIM HD (EYES):

Hard Disk SSD da 128 GB (55.000 immagini) Software Measure incluso DICOM disponibile come opzione.



IMPROVING PERFORMANCES